



# THE PROMISE LIGHT WAND



## The Promise Light Wand Instruction Manual

*Product may not look exactly the same as the one shown in this manual.*



## Unpacking Instructions

Please read this manual thoroughly before operating this unit and retain it for future reference.

Immediately upon receiving the product, carefully unpack and check the contents to ensure all parts are present. Notify the shipper immediately and retain packing material for inspection if any parts appear damaged from shipping or mishandling. Save the package and all the materials. In the event that a product must be returned to the factory, it is important to do so in the original box and packing.

## Warning

**You understand that these products are for educational and research use and not as substitute for medical advise. Follow laser safety procedures of never look directly into the beam or its reflection. Always use laser through crystal spheres only not directly pointed to the skin.**

**INFRARED LASERS can be hazardous and have unique safety considerations. Permanent eye injury and blindness is possible if lasers are used incorrectly. Pay close attention to each safety REMARK and WARNING statement in the manual.**

**This manual contains important laser system safety and operation information. Read and understand all instructions prior to powering on laser unit the first time, to avoid laser eye injury and to avoid breaking the law.**

**This product is in perfect condition before delivery, please check all accessory when unpacking the unit. User should follow this manual to operate the product and make sure that everything is clear on all laser safety and product operation information before powering it on. Any wrong operating will make this unit out of warranty.**



## Packing list

The following items are supplied:

1. Promise Light Wand with Golden Pyramid head
2. Tripod
3. Protective glasses
4. Small Crystal Sphere Set
  - a. Clear Quartz – Photoporation
  - b. Amethyst- Crown Chakra
  - c. Lapis Lazuli - Brow Chakra
  - d. Aquamarine - Throat Chakra
  - e. Green Aventurine - Heart Chakra
  - f. Yellow Calcite - Solar Plexus Chakra
  - g. Carnelian - Sacral Chakra
  - h. Black Tourmaline - Base Chakra



Seven beads resonating with the 7 chakras and their associated organs.



## Operating Instructions

### **Important Note:**

**The Promise Light Wand uses Infrared lasers which are almost invisible to the eye. Never operate The Promise Light Wand without the provided infrared glasses.**

**Never look at the beam directly or to its reflection with or without provided glasses.**

The Promise Light Wand power switch is on the silver end of the wand, it is Off when the switch is pushed in and it is On with a slight push, *switch should appear longer in the On position.*


You can see a slight dim red light when the Promise Light Wand is on by pointing it on your palm.

**NEVER LOOK DIRECTLY INTO THE LASER APARTURE TO VERIFY THE OPERATION.**




## Amplified Crystal Energy

1. Remove the golden pyramid head by twisting away the magnetic locks and pulling out from the wand.



2. Place a small Crystal Sphere of choice into the pyramid head and align its magnetic locks with the Promise Light Wand.



3. When all is assembled start the projection pattern and expose to the area of treatment.





## Photoporation Sphere Amplified Remedy Energy

1. Follow the same steps as shown above; use a dropper or syringe to place several drops of a substance such as remedy, flower essence, herbal tincture etc. into the Clear Photoporation Crystal and seal with silicon caps tightly.



2. Place remedy filled Clear Photoporation Crystal into the freezer.
3. After freezing, place remedy filled Clear Photoporation Crystal into the golden pyramid head and press the power button and expose to area being treated.



## Trouble-shooting instructions

[www.ThePromiseRevealed.com](http://www.ThePromiseRevealed.com)



Faults	Causes	Solutions
Cannot turn on	Battery is low.	Replace the batteries.
	The end cap isn't tightened.	Tighten the end cap.
	Batteries inserted incorrectly	Battery + side facing up towards the end cap

## Notes and Contraindications

*Use only with Crystals.*

*Never look directly at laser beam or its reflections.*

*Never direct beam without Crystal to the skin.*

### Do not use on:

1. Pregnant women
2. Cancer
3. Hemorrhagic diseases

## Maintenance & Cleaning

1. Use cotton soaked in neutral solution to clean the watch surface. Do not use organic solution or gasoline.
2. Please keep away from humidity, heat source and direct exposure to sunlight, in order to avoid the degradation of components which may shorten the service life of the watch.
3. Please do not place the wand with hard objects in order to avoid scratches on its surface or damage of the laser.
5. Please use good quality batteries to avoid leakage that can damage the Promise Light Wand.

## Product Warranty

The Promise Light Wand, will provide free service warranty for one year against manufacture defects.

Our company will charge repair fees for the followings excluded in the scope of repair guarantees;

1. Damages caused by improper operation or mistakes.

[www.ThePromiseRevealed.com](http://www.ThePromiseRevealed.com)



2. Breakdowns caused by self-disassembling.
3. Material costs should be charged for products over the warranty period.
4. Breakdowns or damages caused by accidents or human activities (dropping, scratches, soaking or self-replacement of parts).

## Technical Specifications

Laser medium: GaAlAs semiconductor diodes

Laser wavelength: 808 nm

Maximum power output of single terminal laser: 500 mW

battery voltage: 2 X AAA 1.5V

Ambient temperature: 5-40 °C

Relative humidity: ≤ 80%